

Make Joint Prep Simple!

— THE —[®] CARTILAGINATOR

A revolutionary joint preparation technology designed to precisely prepare the indicated joint in a fast and more efficient manner

X PRECISION

Engineered to control the depth of joint preparation and **minimize bone removal**

X EFFICIENCY


Designed to **decrease joint prep time** compared to traditional joint prep tools*

X CONFORMITY


Anatomically contoured to match the joint, facilitating congruency of the distal and proximal joint surfaces

X EASE OF USE

Allows for a **minimal incision** surgical approach



*Large
Concave*



EXTREMITY[®]
MEDICAL

Flat
Single/Double Sided



Flat
Saw Tip



Large

Curved (3D)
Concave/Convex



Small

*Procedure Type	Flat Single-Sided	Flat Double-Sided	Large Curved, Concave (3D)	Large Curved, Convex (3D)	Small Curved, Concave (3D)	Small Curved, Convex (3D)	Flat Saw Tip
	168-05010	168-05020	170-01250	170-01251	170-02140	170-02141	170-00011
Ankle Fusion	✓	✓	✓	✓			✓
Charcot	✓	✓					✓
Double Arthrodesis	✓		✓	✓	✓	✓	
Lapidus	✓	✓					✓
Midfoot Fusion	✓	✓					✓
MTP Fusion	✓				✓	✓	
NC Fusion	✓	✓	✓	✓			
Subtalar Fusion	✓	✓	✓	✓			✓
Triple Arthrodesis	✓				✓	✓	✓

Part Numbers include “-S” and Descriptions include “,Sterile” for Sterile Packaged

*Based on surgeon reported usage

973.588.8980
ExtremityMedical.com | TheCartilaginator.com
customerservice@ExtremityMedical.com

300 Interpace Parkway, Building A, Floor 2 | Parsippany, NJ 07054

Extremity Medical® and The Cartilaginator™ are trademarks of Extremity Medical, LLC.
© 2024 Extremity Medical, LLC. All Rights Reserved.

LBL-170-00003-EN REV A

LAVENDER MEDICAL

Lavender Medical Ltd.

11 Amor Way, Letchworth Garden City,
SG6 1UG, UK
Tel: 0345 676 9733
Email: info@lavendermedical.com
www.lavendermedical.com

