


SILASTIC TENDON IMPLANTS

FLEXOR TENDON RECONSTRUCTION



Distributed by Lavender Medical, these are
AVAILABLE NOW IN THE UK



The temporary silicone tendon replacement is a monolithic die-cast body with the advantage of an oval cross-section and elongated shape, that comes in a range of diameters and lengths to suit.


ORDERING CODE	DESCRIPTION (W X H X L)
Z06230-0000-03120	Temporary silicone tendon replacement 3mm x 1.5mm x 120mm
Z06230-0000-03180	Temporary silicone tendon replacement 3mm x 1.5mm x 180mm
Z06230-0000-03220	Temporary silicone tendon replacement 3mm x 1.5mm x 220mm
Z06230-0000-04120	Temporary silicone tendon replacement 4mm x 2mm x 120mm
Z06230-0000-04180	Temporary silicone tendon replacement 4mm x 2mm x 180mm
Z06230-0000-04220	Temporary silicone tendon replacement 4mm x 2mm x 220mm
Z06230-0000-05120	Temporary silicone tendon replacement 5mm x 2.5mm x 120mm
Z06230-0000-05180	Temporary silicone tendon replacement 5mm x 2.5mm x 180mm
Z06230-0000-05220	Temporary silicone tendon replacement 5mm x 2.5mm x 220mm
Z06230-0000-06120	Temporary silicone tendon replacement 6mm x 3mm x 120mm
Z06230-0000-06180	Temporary silicone tendon replacement 6mm x 3mm x 180mm
Z06230-0000-06220	Temporary silicone tendon replacement 6mm x 3mm x 220mm

Orders can be placed by contacting orders@lavendermedical.com

Distributed by Lavender Medical, these are
ALSO AVAILABLE IN THE UK



SILASTIC MCP/PIP REPLACEMENT



APEX: IP FUSION DEVICE



**XMCP™: INTRAMEDULLARY
FUSION DEVICE**



**NEOSPAN STAPLES
INDICATED FOR THE HAND**